


AD-A124 782


IMPLEMENTATION OF SOFTWARE FOR A MICROCOMPUTER BASED  
MOSSBAUER SPECTROMETER(U) AIR FORCE INST OF TECH  
WRIGHT-PATTERSON AFB OH SCHOOL OF ENGINEERING B E PATE  
MAR 82 AFIT/GNE/PH/82M-11 F/G 9/2

1/1

UNCLASSIFIED

NL





MICROCOPY RESOLUTION TEST CHART  
NATIONAL BUREAU OF STANDARDS-1963-A

AFIT/GNE/PH/82M-11

AD A124182

IMPLEMENTATION OF SOFTWARE  
FOR A MICROCOMPUTER BASED  
MOSSBAUER SPECTROMETER

Thesis

AFIT/GNE/PH/82M-11      Bruce E. Pate  
                                 Captain    USA

DTIC FILE COPY

DTIC  
ELECTE  
FEB 23 1983  
S D E

Approved for public release; distribution unlimited.

IMPLEMENTATION OF SOFTWARE FOR A  
MICROCOMPUTER BASED  
MOSSBAUER SPECTROMETER

THESIS

Presented to the Faculty of the School of Engineering  
of the Air Force Institute of Technology  
Air University  
In Partial Fulfillment of the  
Requirement for the Degree of  
Master of Science

by

Bruce E. Pate, B.S.  
Captain USA

Graduate Nuclear Engineering

March 1982

Accession For	
NTIS GRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/ _____	
Availability Codes	
Avail and/or	
Dist	Special
A	



Approved for public release; distribution unlimited.

## Preface

This thesis is the product of my efforts to program an PDP-11 microcomputer to analyze the spectra generated by a Mossbauer spectrometer, so that a more responsive, interactive, and versatile evaluation can be made. This thesis will describe the program used to evaluate the spectra, the modifications made to the program to better work with the current spectrometer, the results of test evaluations made using absorbers of potassium ferrocyanide, natural iron, sodium nitroprusside, and neodymium cobalt, and how these results compared with the previous computer program and with published standard values.

I would like to thank Dan Zambon of the Electrical Engineering Department, whose assistance in the initial set up of the computer was invaluable. I especially want to thank my advisor, Dr. George John, whose guiding hand and support during this project helped make this thesis possible.

I also want to thank my wife, Kathy, whose moral support and understanding were always present when I needed them.

Bruce E. Pate

## Contents

	<u>Page</u>
Preface . . . . .	ii
List of Figures . . . . .	v
Abstract . . . . .	vi
I. Introduction. . . . .	1
Problem . . . . .	1
Background. . . . .	2
Theory. . . . .	3
Scope . . . . .	5
II. Equipment Used. . . . .	7
Mossbauer Equipment . . . . .	7
Source and Absorbers. . . . .	8
Computer System . . . . .	9
III. Data Processing . . . . .	10
Mossbauer Parameters. . . . .	10
Program GENFIT. . . . .	13
MOSFUN Background . . . . .	14
Main Program. . . . .	15
Subroutine Descriptions . . . . .	17
DATI. . . . .	17
DRIVE . . . . .	18
PARM. . . . .	19
TMIN. . . . .	19
SAVE. . . . .	20
NORM. . . . .	20
FIT . . . . .	20
FTHEO . . . . .	22
THEO1 . . . . .	22
THEO2 . . . . .	23
SUBT. . . . .	24
SHIFT . . . . .	24
TRANS . . . . .	24
LINEQ . . . . .	25
PRNT. . . . .	25
PLOT. . . . .	25
PUNC. . . . .	25
Analysis Commands . . . . .	26
IV. Results . . . . .	30
General Procedure . . . . .	30
Potassium Ferrocyanide. . . . .	32
Natural Iron . . . . .	33
Sodium Nitroprusside. . . . .	33

	<u>Page</u>
Neodymium Cobalt. . . . .	34
Program Considerations. . . . .	35
V. Conclusions and Recommendations . . . . .	40
Bibliography. . . . .	41
APPENDIX A: MOSFUN Installation Procedure . . . . .	42
APPENDIX B: Subroutine DATI . . . . .	45
APPENDIX C: Subroutine DRIVE. . . . .	50
APPENDIX D: Subroutine PLOT . . . . .	56
APPENDIX E: MOSFUN User Supplement. . . . .	59
Vita. . . . .	73

### List of Figures

<u>Figure</u>		<u>Page</u>
1	Example of a Mossbauer Spectrum . . . . .	4
2	Natural Iron Results for Peaks 3 and 4 - MOSFUN Using Lorentzian (Areas). . . . .	37
3	Sodium Nitroprusside Results - MOSFUN Using Voigt (Inten.). . . . .	39



Abstract

Computer analysis of Mossbauer spectra using a micro-computer requires tight controls on memory requirements and programming techniques. In this study, a FORTRAN IV program, MOSFUN, was programmed into an PDP-11 microcomputer. The program was modified to calculate the velocity for a data point using the least squares fit of counts from a Moire interferometer and to perform an alphanumeric plot of the original data and the theoretical fitted spectrum. Spectra were obtained for absorbers of potassium ferrocyanide, natural iron, sodium nitroprusside, and neodymium cobalt. These spectra were then analyzed with the MOSFUN program and the results compared with analyses done with the former program, GENFIT, and with published values. Results showed MOSFUN to be more flexible, easier to adapt, and easier to control during the evaluation of spectra, while still providing results that agreed within 1% of those obtained with GENFIT and with standard values.

# IMPLEMENTATION OF SOFTWARE FOR A MICROCOMPUTER BASED MOSSBAUER SPECTROMETER

## I. Introduction

### Problem

The ultimate purpose of this project is the design, construction, and testing of a microcomputer-based system capable of collecting, storing, and analyzing data from a Mossbauer spectrometer. This new system will replace an aging and unreliable multichannel analyzer (MCA) based system and will preclude the requirement to do analysis on a separate mainframe computer. This new system will allow data to be collected more efficiently and will provide more productive and timely analysis.

The objectives of this thesis were to program an PDP 11 microcomputer to analyze Mossbauer spectra and to assist in the design of new hardware to replace the current MCA. The software will allow more control over the fit of the Mossbauer spectra, since it will be interactive instead of the batch process of the current program. Analyses can be done more quickly and frequently since the computer is dedicated to the spectrometer. The new hardware is designed to provide the signals needed by the Mossbauer control unit to drive the velocity transducer, to multiplex the incoming signals, and to store the data until requested by the computer.

## Background

Until recently, Mossbauer spectrometers have been based on systems which used multichannel analyzers (MCA) for acquiring and storing spectra. These systems required the data to be collected and then transferred in some way to a computer for analysis of the data. The transfer can be accomplished in a variety of ways, which vary in speed and tedium. The most rapid method would be a direct transfer of the data stored in the MCA memory to the memory of a computer. This may be by a direct connection to a minicomputer or over phone lines to a large centrally-based time-shared computer. An intermediate, but more time-consuming method, would be to transfer the data to magnetic tape or paper tape before transferring to a computer. In the case of paper tape and computer systems with inadequate input directly from paper tape, another step requiring the transfer of the data to punched cards may be required. The later example with paper tape has been the situation with the Mossbauer system and the facilities for analysis available at this laboratory.

An attractive alternative to the foregoing is the use of a microcomputer directly interfaced to a Mossbauer transducer and driver. This alternative has been enhanced by hardware and software which in recent years have become available in large variety and at moderate cost. In addition, programs for analysis of Mossbauer spectra have become available for a nominal fee. One such program, called MOSFUN, has been written in FORTRAN IV by E.W. Muller (Ref 3:1.1).

## Theory

Mossbauer spectroscopy makes use of the Mossbauer effect, which was discovered by Rudolf Mossbauer in 1957 (Ref 2:12). By moving a radioactive source, Mossbauer discovered that the energy of the emitted photon was changed by a small amount. The energy was increased for a positive velocity and decreased for a negative one. This effect is used for spectroscopy by placing an absorber between the moving source and a detector. The absorber to be studied must contain the same atom as the source but in the ground state. If the energy of a photon exactly matches the energy needed to raise a nucleus in the absorber to a higher state, the photon may be absorbed. The number of counts seen by a detector will be reduced whenever the energy of a photon emitted by the source is changed enough to equal this resonance absorption energy. A plot of the number of counts recorded by the detector versus the velocity of the source results in a spectrum such as the one in Figure 1. The dips in the spectrum correspond to the excited energy levels in the absorber nuclei. The velocity of the source can be related to the energy of the absorbed photon and, from these energies, a great deal can be told about the states of the nuclei in the absorber. Each dip in the spectrum may be one or more overlapping energies. Each energy is called a line and each line will have a common shape. The shape is a convolution of the energy distributions from the source and the absorber. This distribution will be




Figure 1. Example of a Mossbauer Spectrum

Lorentzian unless perturbed by the system. Detailed descriptions of Mossbauer theory and its application to spectroscopy can be found in the thesis by Towery (Ref 11:6).

A computer program to analyze this spectrum must be able to apply a Lorentzian distribution to each line, determine the intensity or magnitude of each dip, and locate the velocity of the lines using several known characteristics or parameters of the nuclear energy levels. The program then varies the parameters until it finds a set of theoretical lines, whose combined shapes and locations most nearly match the experimental data.

### Scope

To interface the microcomputer to the Mossbauer system requires interface hardware to be built and tested to replace the functions currently done by the MCA. This new device must be interfaced with the PDP 11 microcomputer and the computer programmed to use the interface as an extension of its own memory. The new device has not been completed and these parts of the project will not be covered in this report.

Once in the computer, the data must be analyzed. This study is primarily concerned with the computer programming needed to analyze Mossbauer spectra after the data is in the computer. Chapter II discusses the equipment used in the spectrometer and the computer system used to evaluate the spectra produced. Chapter III discusses the parameters used to evaluate Mossbauer spectra, the computer program

MOSFUN and how it uses the parameters to evaluate spectra. Chapter IV discusses the results of experiments done with natural iron, sodium nitroprusside, potassium ferrocyanide, and neodymium cobalt, and how these results obtained using MOSFUN compared with results obtained with the current program, GENFIT, and with published standards. Chapter V discusses the conclusions and recommendations for this study.

## II. Equipment Used

### Mossbauer Equipment

In addition to the source and absorber, the components of the Mossbauer spectrometer are a velocity transducer (motor), a motor control unit (MCU) to drive the motor at a constant acceleration, a linear amplifier/single channel analyzer, a krypton-filled proportional counter, and a high voltage power supply for the detector. All except the detector were manufactured by Ranger Electronics. The MCU is driven by the channel-address ramp from an RIDL multichannel analyzer (Model 34-12B) driven in the time-sequential scaling mode by an RIDL time-base oscillator (Model 54-6). The equipment is the same as that used by Beller (Ref 1:10), and described in detail in the thesis by Skluzacek (Ref 8:5-9), and will not be repeated here.

The transducer can drive the source in either of two drive modes, triangular or flyback. The maximum velocity (VMAX) of the source is set on the MCU and the source is moved at a constant linear acceleration between +VMAX and -VMAX. The triangular drive mode is a full period drive mode and, as such, moves the source from +VMAX to -VMAX and back to +VMAX. It continues to cycle the source through this full period as counts from the detector accumulate in the MCA. Each pass between +VMAX and -VMAX



generates a complete spectrum such as the one in Figure 1. As a result, the full period drive mode generates two spectra, one from +VMAX to -VMAX and a mirror image spectrum from -VMAX to +VMAX. The flyback mode on the other hand is a half period drive mode. The source velocity is only varied from -VMAX to +VMAX and generates only a single spectrum rather than the double spectrum of a full period drive mode. The spectrum in Figure 1 is a half period spectrum (flyback) and is identical to the right half of a full period spectrum such as a triangular drive mode would create. The flyback drive mode was used exclusively in this study, but references will be made to programming considerations for full period spectra.

#### Source and Absorbers

The source used for obtaining spectra used in this study is 50 mCi cobalt-57 (date: 5 June 1981) diffused into a 6 micron thick rhodium matrix. The source is prepared by Amersham-Searle. The cobalt-57 decays into an excited state of iron-57, which in turn decays emitting an 14.4 kilo-electronvolt (keV) gamma ray among other radiations. This 14.4 keV gamma ray is the photon of interest for this study. Therefore, the absorber used must also contain iron-57 in the ground state that can be excited with a photon of this energy. Absorbers used in this study are a 25 micron thick natural iron foil, a National Bureau of Standards sodium nitroprusside crystal, a potassium ferrocyanide crystal, and a neodymium cobalt sample.

### Computer System

The computer system consists of a PDP 11/03 micro-computer, an RX-02 dual 8-inch floppy disk drive, and an VT-100 terminal all manufactured by Digital Equipment Corporation. The memory is a 64 kilobyte CI-1103 memory board made by Chrislan Industries. An MX-80F/T dot matrix printer by Epson is connected serially to provide printouts of results.

The PDP 11/03 is a 16 bit word computer. It consists of an LSI 11/02 microprocessor with an KEV-11 floating point chip, an BDV-11 bootstrap and diagnostic board, an DLV-11 4-port serial interface board, and an DRV-11 parallel interface board. The operating system is RT-11 version 04B. Programming is done in FORTRAN IV version 2.5. The FORTRAN software is also provided by Digital Equipment Corporation.

### III. Data Processing

The analysis of Mossbauer spectra makes use of several characteristics of the Mossbauer effect. These characteristics are called parameters and the parameters have the same relationships regardless of the absorber used. These parameters will determine the relative intensity of each line to the other lines, the half-width at half-maximum of the dip corresponding to a line, and the velocity position of the line relative to zero velocity. A computer program to analyze the spectra must be able to use these parameters to compute a theoretical spectrum that most nearly matches the experimental data.

#### Mossbauer Parameters

The parameters that are used to compute the theoretical spectrum (or fitting the spectrum) are divided into three groups: intensity, position, and half-width. If the energy distribution is not truly a Lorentzian, a fourth set of parameters called form may be required.

The intensity of a line is defined as the number of counts between the lowest portion of a dip and the baseline. The baseline is the horizontal line across the top of the spectrum that represents the spectrum from the source photons if none were absorbed (see Figure 1). For example, if the baseline contained 100000 counts and the lowest part of a dip contained 80000 counts, the intensity of that line would

be 20000 counts. This requires that the baseline and the intensity of each line be variable parameters. Intensities of some lines are related so that one intensity parameter may be used for more than one line. The intensity of a line may also be expressed as an area, which is the intensity multiplied by the half-width at half-maximum. In some absorbers (natural iron for example), the areas of two lines may be closer related than the intensity and the area may be the best parameter to use.

Since the intensity depends strongly on accurately determining the baseline, some additional parameters may be required. In reality, the baseline is not flat, but has a slight curvature caused by the varying source-to-absorber distance. The closer the source is to the absorber, the higher the number of counts in the baseline. This difference is caused by the increased solid angle the absorber represents when the source is nearer. This curvature may be handled in several ways including polynomial fits to the baseline, allowing the geometry to be a variable, or assuming the baseline is a straight line. If the spectrum to be fit is a full period spectrum, there is a point of symmetry in the spectrum. This symmetric point, called the fold point, represents the point of inflection in the baseline curvature. The baseline is concave in one half of the spectrum and convex in the other. If the spectrum is folded in half at this point, the baseline curvature will tend to cancel out and the dips in each of the mirror spectra will

coincide, increasing the statistics of each dip. The curvature and the fold point must be considered when determining parameters to use to fit Mossbauer spectra.

The position of a line relative to zero velocity (and thus relative to the unshifted source photon energy) is determined by some combination of three parameters; isomer shift, hyperfine magnetic splitting, and quadrupole coupling. The isomer shift is the amount the midpoint of the spectrum is shifted from zero velocity. The isomer shift is caused by the difference in the nuclear radii of the ground and excited states, and how this difference couples with the electron charge density at the nucleus. The hyperfine magnetic splitting results from a coupling of the nuclear magnetic moment with the magnetic field due to the atom's own as well as adjacent atom's electrons. The result will be pairs of lines corresponding to transitions between nuclear magnetic sub-levels with the lines in each pair separated by an amount proportional to the magnetic field present. Quadrupole coupling results from the interaction of the nuclear quadrupole moment with the gradient of the electric field produced by other charges in the crystal. The quadrupole coupling also results in a splitting of lines. These effects and the Hamiltonians used to represent them are described in detail by Roberts (Ref 5:5). The position of a line may be determined by a combination of these effects but the combination is consistent for a particular type of absorber.

The half-width parameter is the half-width of a dip at

half the maximum intensity of the dip. It is measured in velocity units and is related to the Lorentzian lineshape of the dip.

The form parameters are only used if the distribution function to be used is not pure Lorentzian. These parameters may be used to change the lineshape from Lorentzian to Gaussian or to some combination of both to account for fluctuations in the energy distribution of a line.

#### Program GENFIT

The computer program GENFIT is a least-squares minimization curve fitting program developed at Argonne National Laboratory. This has been the computer program used in the past to analyze Mossbauer spectra. The design and memory requirements of GENFIT require a mainframe computer and, therefore, this program could not be used on the PDP 11 micro-computer. It is, however, used as a reference to compare with the results obtained using MOSFUN and the PDP 11 system.

GENFIT uses a user-supplied subroutine called CALFUN to relate the Mossbauer parameters to the spectrum from a particular absorber. This subroutine is the model used by GENFIT to find the best fit to the experimental data. GENFIT assumes a quadratic curvature to the baseline and then subtracts this curvature as background from the baseline. This leaves a linear baseline for intensity calculations. The shape, relative intensities, position, and half-width of the lines is determined by the model described in the CALFUN

subroutine to be written for each absorber type used. In addition, no provisions exist to modify the model during an evaluation. GENFIT is described in more detail by Skluzacek (Ref 8:16).

#### MOSFUN Background

The computer program used for this study is a FORTRAN IV program called MOSFUN. It was written by E.W. Muller at the Johannes Gutenberg-Universitat in Mainz, West Germany and is available through the Mossbauer Effect Data Center, University of North Carolina (Ref 4:90).

MOSFUN is written in a very general manner so that it can be used on a variety of computers. It is available in three formats: batch, console (interactive), and a modified console version for use on PDP systems. Since the program is written very generally, there are no machine-specific functions used and no special spectrometer capabilities required.

Four models for use in fitting spectra are provided, with simple provisions for expanding to more complex models. Two models (called theories by MOSFUN) are done using Lorentzian lineshapes, one using dip intensity and one using dip area. The other two theories use approximate Voigt lineshapes, again one using dip intensity and one using dip area. Iterations can be carried out using either Newton's method or the Gradient method to determine the corrections required in the theoretical spectrum.

MOSFUN uses the geometry of the baseline to account

for its curvature and is capable of analyzing both full period and half period spectra. The parameters are read from a disk file rather than requiring new subroutines to change absorbers. The parameters can be modified during the fit process if one of the interactive versions is used.

The PDP version of MOSFUN is used for this study with some modifications to better use the current spectrometer. The PDP version is basically the console version (designed to run on mainframe computers using terminals) reduced in size, so that it can be overlaid into the PDP computer's smaller memory. MOSFUN makes extensive use of subroutines and this makes overlaying much easier. MOSFUN contains 17 subroutines that are divided into two overlay regions. The PDP version allows up to 32 fit parameters, 24 lines, and 1024 data points (Ref 3:4.14). Three of the subroutines (DATI, DRIVE and PLOT) were extensively modified, as was the main program in the areas that used these subroutines. Installation of MOSFUN on the PDP 11/03 is described in Appendix A.

#### Main Program

The main program for MOSFUN (Ref 3:3.3) acts primarily as an interactive guide. One of 16 commands (Ref 3:4.16) is entered from the console, then the main program determines what initial values are needed and which subroutines are to be called. Initial values are requested and entered from the terminal. Data and parameter values are stored on disk files and are entered from the disk. Variables are transferred between the subroutines and the main program through common



blocks. Output can be directed to the terminal or the line-printer.

The PDP 11 uses a technique called overlaying to allow programs much larger than the computer's memory to execute. In this process, more than one subroutine can share the same computer memory space as long as both are not needed at the same time. MOSFUN takes advantage of this capability by placing all subroutines that execute a command in one overlay region. Any subroutine that these subroutines may need to call is placed in a second overlay region. If more than one subroutine is needed at the same time, one may be placed in the main root segment which is always in the computer's memory. A breakdown of the overlay regions and the subroutines that each contains is:

<u>Region</u>	<u>Subroutines</u>				
Main Root	MOSFUN	DRIVE	FTHEO		
Overlay Region 1	DATI	PARM	TMIN	SAVE	
	FIT	NORM	SHIFT	SUBT	
	PRNT	PUNC	PLOT		
Overlay Region 2	TRANS	LINEQ	THEO1	THEO2	

When a command is given, the main program, MOSFUN, calls the appropriate subroutine from overlay region 1. These subroutines handle the experimental data, the parameters, the correlation matrix, the fit procedure, spectrum manipulation, and the output of results. The subroutines in overlay region 2 enter any values into the program from the console, solve a

system of linear equations, and model the theoretical spectrum. Each command is entered interactively from the terminal by typing the first two letters of the command. The commands may be given in any order after the experimental data is entered and the parameter starting values described.

### Subroutine Descriptions

A brief description will be given concerning the function of each subroutine. Each subroutine and its variables and indicators are described in detail in the MOSFUN documentation (Ref 3:3.1-3.27), so will only be reviewed here. The modifications made to the program will be covered in more detail. The word "channel" describes the address of a data point within a series of storage locations. This may be in an MCA or a computer memory chip.

Subroutine DATI (Ref 3:3.8) handles the input of experimental data from the disk and sets the channel limits over which the spectrum will be fit. The experimental data is checked point by point for dropped counts and DATI provides the opportunity to zero any channel that appears bad. A bad data point is one that differs from the previous channel by more than statistical error times a correlation factor entered from the console. Any channel that has been zeroed is ignored for fit calculations and is not plotted. Any channel outside the fit limits is set negative and is ignored for fit calculations but is plotted. Overflows are added to the data during this subroutine.

DATI was modified to identify and sort out counts

from the internal time-keeping oscillator and the Moire interferometer for velocity measurement which are multiplexed with the data by the spectrometer. The counts, proportional to velocity when extracted, are corrected for overflows, converted to a velocity, and related to a channel number by a linear least squares fit. These data points are then set to zero so that the remaining data are identified as spectral data. Appendix B contains a listing of the changes to DATI.

Subroutine DRIVE (Ref 3:3.22-3.26.3) calculates the velocity and geometry factor for a particular channel. It can also calculate the derivatives for parameters 1 (baseline), 2 (geometry), and 3 (folding point) which are used during the fit process. DRIVE calculations depend on the drive mode used with the transducer. Triangular wave and sine wave modes are provided in MOSFUN. The geometry factor accounts for the curvature of baseline caused by the varying source-to-absorber distance. The fold point represents the point of symmetry in a full period spectrum such as that generated by a triangular wave or sine wave. The velocity in any channel is determined by the maximum velocity (VMAX) relative to the channel's position in the period.

DRIVE was modified to calibrate the velocity relative to a channel from the velocity data extracted during DATI. This velocity fit is used to calculate the velocity for any channel. A third drive mode was added since all current work is done using the flyback mode. The flyback spectrum is essentially the right half of a triangular wave spectrum and

computations for geometry were carried out as such. The fold point is not used for flyback spectra since flyback is only a half period spectrum. Even though the fold point is a parameter it does not enter into the calculation of half period theoretical spectra. Appendix C contains a listing of the changes to DRIVE.

Subroutine PARM (Ref 3:3.12-3.13) controls the input of, and changes to, the fit parameters. Each parameter is defined by an index number, a name, a value, and upper and lower limits. Once defined, the number is the only quantity required to reference that parameter. The upper and lower limits are bounds which restrict the values the parameter can assume. If the parameter exceeds either limit, the parameter is set equal to that limit for the next iteration and a warning message is displayed. If the limit is set equal to zero, the program is free to adjust the parameter to any value. Each parameter can be set to vary (FREE) or can be treated as a constant (FIX). Any time a command can be given, a parameter can be fixed or freed. This allows control over well defined parameters during the first few blocks of iterations. Four values of the parameters are available. The current values of these parameters can be changed to the starting values, the optimum values (the lowest Chi-squared), or the values of the previous iteration. The parameters must be replaced as a set. This would be done if divergence occurs or if convergence to a local minimum is occurring.

Subroutine TMIN (Ref 3:3.9-3.10) controls the correla-

tion matrix. This matrix is not a correlation matrix in the usual sense of error correlations, but rather correlates the lines in a spectrum to the parameters being used to fit that spectrum. This matrix determines how much of each parameter is used to calculate each line. It is actually four sub-matrices, one each for intensity, position, half-width, and form. Each submatrix has the dimensions of the total number of lines by the number of parameters that affect that sub-matrix. For example, if there are four lines in a spectrum and two intensity parameters, then the intensity submatrix will be a four by two matrix. The values in the correlation matrix are constants and do not vary during the fit of a spectrum.

Subroutine SAVE (Ref 3:3.19) enables the values of the parameters and the correlation matrix to be saved in one disk file. After any change is made to either set of values, the option to change the file (called the "save file") is offered. If changes are saved, the same save file can be used or another file created to save the values, if both sets are desired.

Subroutine NORM (Ref 3:3.14) performs a two parameter fit for the baseline and a scaling factor for the intensity. Relative intensities can then be used during the fit process to calculate real intensities and the baseline. Norm does not calculate Chi-squared or any theoretical information directly.

Subroutine FIT (Ref 3:3.16-3.17) initiates calculation

of Chi-squared and the theoretical spectrum. It performs iterations by either Newton's method or the Gradient method to calculate correction factors to be added to the parameters for the next test cycle.

Chi-squared is used as the test for convergence. It is normalized to the number of fitted channels and the number of free parameters. It is defined as (Ref 3:3.16):

$$\text{Chi-squared} = \left( \sum^N (Y(i) - F(i))^2 / Y(i) \right) / (N - M) \quad (1)$$

where

N is the number of fitted channels

Y(i) is the i-th experimental data point

F(i) is the i-th theory value

M is the number of free parameters

When successive Chi-squared values differ by less than a specified convergence criteria (0.001 by default), iteration ceases and convergence is declared. If a Chi-squared value is larger than the previous value by more than ten times the convergence criteria, the theory is declared to diverge. Iteration stops regardless of the convergence status if the number of iterations specified by the user have been reached.

Newton's method (Ref 3:2.1-2.2) solves a system of equations (using subroutine LINEQ) created during the calculation of Chi-squared. The solutions represent the corrections to the fit parameters. Newton's method will diverge in the vicinity of flat minima such as those of the baseline. Therefore, care must be taken to use only the regions of the

spectrum that contain dips rather than fitting the entire spectrum.

The Gradient method (Ref 3:2.3-2.5) calculates the values of the corrections from the gradient of the previous values of the parameters. It works more slowly than Newton's method but will usually converge. A path width is specified which limits the amount the parameters will be changed in a single iteration. A large path width allows larger changes to be made in the parameter's value and may allow divergence if care is not taken. Too small a path width may delay convergence. The default value of the path width (1000) was used successfully in this study.

Subroutine FTHEO (Ref 3:3.27) resides in the main root region and directs which theory model will be used to calculate the theoretical spectrum. FTHEO is called by subroutine FIT and in turn calls the appropriate theory-model subroutine from overlay region 2 based on the theory index number selected during parameter entry.

Subroutine THEO1 (Ref 3:5.4-5.6) calculates the parameters using simple correlated Lorentzian curves. If a spectrum is composed of  $M(i)$  Lorentzians of intensity  $W(i)$ , position  $E(i)$ , and half-width  $H(i)$ , then the theoretical value  $F$  at velocity  $v$  is given by (Ref 3:5.4):

$$F(v) = g \left( B - \sum_{i=1}^{M(i)} \left( W(i) \frac{1}{1+a^2} \right) \right) \quad (2)$$

where

$g$  is the geometry factor at velocity  $v$

$B$  is the baseline

$$a = (v - E(i)) / H(i)$$

The correlation matrix entered during TMIN determines how the parameters are combined to determine  $W(i)$ ,  $E(i)$ , and  $H(i)$ .

Subroutine THEO2 (Ref 3:5.8-5.12) calculates the theoretical values using Voigt lineshapes. A Voigt profile is the result if statistical disturbances change the Lorentzian lineshape. A Voigt lineshape is a convolution of a Gaussian and a Lorentzian shape. This function cannot be analytically determined but can be approximated by the biquadrated (fourth power polynomial) function (Ref 3:5.8):

$$\frac{1}{1 + ax^2 + (1-a)x^4} \quad (3)$$

A simple Lorentzian is produced when  $a=1$  and  $a=0.5$  approximates a Gaussian; thus " $a$ " is constrained to be between 1 and 0.5. The theoretical value  $F$  at velocity  $v$  for a spectrum of  $M(i)$  lines with intensity  $W(i)$ , position  $E(i)$ , half-width  $H(i)$ , and form factor  $A(i)$  is (Ref 3:5.11):

$$F(v) = g \left( B - \sum_{i=1}^{M(i)} W(i) \frac{1}{X2} \right) \quad (4)$$

where

$$X2 = 1 + (A(i) X1^2) + ((1-A(i)) X1^4)$$

$$X1 = (v - E(i)) / H(i)$$

$g$  = geometry factor at velocity  $v$



B = baseline

The line parameters W, E, H, and A again are correlated by the correlation matrix as in THEO1 except a submatrix is added for the form factors.

Subroutine SUBT (Ref 3:3.14) subtracts out the intensity contribution from any line by subtracting the theoretical fit for that line from the original data. This routine may be useful to remove strong well-defined peaks that may be hiding other peaks. The revised spectrum may be plotted or punched for examination, however, the original spectrum is restored before any additional fitting or testing is done.

Subroutine SHIFT (Ref 3:3.18) shifts a full period spectrum to correct for phase errors in the drive system. This places the zero velocity point directly in the center of a folded spectrum. For example, if the zero velocity channel of a folded spectrum containing 400 data points is in channel 204, then SHIFT will adjust the entire spectrum to put the zero velocity in channel 200.

Subroutine TRANS (Ref 3:3.20) allows format-free reading of numbers and is called anywhere in the program where numerical input is required. One input record (up to 80 columns) is entered from the terminal and transferred to TRANS via common blocks. This record can contain up to 16 values in any format separated by a blank. Each number is deciphered separately and then transferred back via common blocks in the vector OUT. The variable IOUT counts how many numbers have been analyzed and is transferred back in the subroutine statement.

Subroutine LINEQ (Ref 3:3.27) solves a system of linear equations by the Cholesky method. The Cholesky method works fast and requires no memory space for provisional results. This method is only applicable to symmetrical positive definite matrices. However, that is all that is generated during Newton's method or during the least squares fit of the velocity data extracted by subroutine DATI. LINEQ can also calculate the inverse of a matrix when it is needed, such as to calculate the parameter error correlations.

Subroutine PRNT (Ref 3:3.15) sends the results to the lineprinter. Options that may be printed are:

- (1) Original data
- (2) Parameters, parameter errors, and error correlations
- (3) Output from the theory program that gives the intensity, position, half width, and form for each line

Subroutine PLOT is written to plot the experimental data and the theoretical spectrum, if it has been calculated. The plot routine provided in MOSFUN will only work with a VT-55 graphics terminal. Since this system does not have that option, a new alphanumeric plot routine was written. The plot can be directed to the terminal or the lineprinter. Intensities are scaled to the baseline and the plot can be done between any two channels. Appendix D contains a listing of subroutine PLOT.

Subroutine PUNC (Ref 3:3.18) stores the experimental and, if calculated, the theoretical spectra in a disk file for fitting in a later run. The spectra are written to the disk

under a file name requested from the terminal. Subspectra can also be calculated for each intensity parameter and their values are written to the disk separately.

### Analysis Commands

The analysis of a spectrum requires several commands, and the order in which they are executed is somewhat flexible. Data, parameter values, and the correlation matrix must be entered, initial testing of the parameter set performed, spectrum folded or shifted if desired and then the fit procedure begun. Results may be printed, plotted, or saved any time a command can be given. Entering commands requires only the first two letters of the command be typed from the terminal. Commands and their prompts are discussed in detail in chapter 4 of the MOSFUN documentation (Ref 3:4.15-4.32).

Entering data is accomplished by the command DATIN. The storage file for the data and a correlation factor are requested. A correlation factor of 100 is sufficient to check for dropped counts in most spectra. The format for data on the disk is explained in Appendix E. The data is loaded into the computer's memory and a velocity calibration done. Any points that appear bad are listed with an option to zero. The option to zero any range of channels and/or set a range of channels that will not be fit is offered next. Any time a command can be given, DATIN may be commanded and the fit limits changed, without changing the experimental data.

The command PARIN allows entering of a parameter set, changing of any parameters value, and saving of the current

values on the disk. Appendix E contains a detailed explanation of how to enter a new set of parameters.

The command CORR enables initial entry of the correlation matrix or changes to the current values. If the correlation matrix is read from a file, it is part of the same file as the parameters and is read during the command PARIN. Details of the correlation matrix entry are described in Appendix E.

The command NORM should be given after entering the data parameters. This enables the program to use relative intensities rather than actual intensities and baseline.

The command TEST will perform an initial Chi-squared calculation on the initial parameters. If the Chi-squared value is too large, the parameters may need to be changed before the fit is begun.

The commands FOLD and SHIFT modify the experimental data to account for a full period spectrum. These commands would be used before the fit process is started. However, all spectra analyzed during this study were only half period spectra and these commands were not needed.

The fit procedure is begun by issuing the command NEWTON or GRADIENT. Both will request the number of iterations desired and the convergence criterion. A default convergence criterion of 0.001 is used if none is specified. For the GRADIENT command a path width is also requested. The default value of 1000 was used during this study. Larger values may allow faster convergence since larger changes can be made to the parameters, however, larger values may also lead to divergence. Small values

may needlessly delay convergence but may be needed if the default value leads to divergence. The number of iterations is kept small initially and the parameters evaluated after each iteration. Changes in the parameters can be made if problems begin to appear during the iteration process.

The commands RESTORE, FIX, and FREE allow control over the parameters. RESTORE changes the current parameter set to either the starting set, the optimum set, or the set from the previous iteration. These changes are made to the entire set of parameters. Selective changes must be made using the command PARIN. Commands FIX and FREE allow any combination of parameters to be fixed as constants or allowed to vary. Fixing well-defined parameters early insures major adjustments are only made to weaker parameters. Freeing only well-defined parameters allows the fit to be made to these parameters, and then the weaker parameters can be selectively released.

After each iteration, the following values are displayed on the terminal: the current value and previous value of each parameter, the Chi-squared value for the current set of parameters, the difference between this Chi-squared and the previous one, and the percent of the theory points that are within a statistical error's range of the experimental value. Statistical error is defined as the square root of the experimental value.

The commands PRINT and PLOT send the results to the lineprinter. PRINT can print the original data, the parameters alone or with their respective errors and error corre-

lations, and each line's value for intensity, position, half-width, and form from the theory program. PRINT also asks for text of up to 60 characters that can be used to describe the conditions that led to the results. PLOT can generate an alphanumeric plot and send it to either the terminal or the lineprinter. The plot is made between two specified channels and includes original data and the theoretical spectrum, if it has been calculated. The velocity scale is printed on the X-axis and relative intensity on the Y-axis.

The command PUNCH initiates transmission of the theory spectrum or subspectra for each intensity parameter to an output file on a disk. This file can then be used for fitting on a later run.

The command SUBTRACT removes the contribution from any theoretical intensity parameter from the experimental spectrum. This may reveal hidden peaks that are buried under much more intense peaks. After the subtraction, the resulting spectrum can be plotted or punched. Before any other command is executed, the original spectrum is restored by the program.

The command EXIT is the last command given, as it terminates program execution. A new analysis can be started without exiting by giving the commands DATIN and PARIN. Any files that are open are closed before execution is terminated.

#### IV. Results

To test the capabilities of MOSFUN to analyze spectra, four absorbers were chosen as test samples: potassium ferrocyanide, natural iron, sodium nitroprusside, and neodymium cobalt. Potassium ferrocyanide presents a simple, one-line spectrum with the location of the dip determined only by the isomer shift. Natural iron presents a six-line spectrum with the isomer shift and hyperfine splitting determining the position of the lines. Sodium nitroprusside has a two-line spectrum with the isomer shift and quadrupole coupling determining the position of the lines. Neodymium cobalt has a complex spectrum with several overlapping six-line iron spectra. All three position parameters are involved in determining the position of each set of lines and the parameters will be different for each set. Spectra for the first three absorbers were taken during the summer of 1981 using the current MCA based Mossbauer spectrometer. The neodymium cobalt spectrum was taken in 1976 by John (Ref 9:117) but used basically the same system. The data was manually stored on a computer disk where it could be recalled for analysis. The data for sodium nitroprusside and natural iron were also transferred to punched cards and analyzed by GENFIT on the mainframe computer for comparison with MOSFUN.

#### General Procedure

A standard procedure was used to analyze each set of data. The data was obtained and put into the computer. A set

of parameters was chosen based on known quantities for each sample. The data was then analyzed. The results obtained were compared to the results obtained with GENFIT and to published standards.

All analysis runs made during this study were half period spectra using the flyback drive mode. Each set of data was stored in 400 channels. Time-base data was stored in the first 14 channels and velocity counts multiplexed in every channel ending in a 9 (19,29,...,399). Mossbauer data was stored in all remaining channels.

After the data and parameter set were entered, commands TEST and NORM were executed. NEWTON was always tried, but always diverged before a final result could be found. Therefore, results shown were obtained using the Gradient method with the default convergence criterion and path width. Five iterations were tried initially and the parameters examined after each iteration. Then blocks of ten iterations were executed until convergence was reached. The number of iterations and the time required for each iteration depended on the number of Mossbauer lines and the number of free parameters.

Both Lorentzian and Voigt lineshapes were tried with each set of data. The final results reached were essentially the same, but convergence always took much longer using the Voigt lineshape. The similarity in the Lorentzian and Voigt results shows that there is very little statistical fluctuation in the experimental data.

The fit could be made over less than the full 400



channels. The range over which the data is to be fit is called the "fit range". In some cases, the fit range was varied to test the influence this had on convergence. No significant differences were noted in the final results; although convergence was generally faster if the fit was performed only over the range of channels containing the dips.

The results give an error limit for some parameters. The limits are somewhat deceiving between MOSFUN and GENFIT. MOSFUN and GENFIT both compute an error correlation matrix that tracks the propagation of system errors. GENFIT estimates the value of the standard deviation of a parameter as the square root of the diagonal element of this error matrix that corresponds to that parameter. The GENFIT error is for 67% confidence so is just this standard deviation. MOSFUN estimates the standard deviation as the square root of the diagonal element multiplied by the Chi-squared value. MOSFUN computes errors on a 95% confidence level which is 1.96 times the standard deviation. Therefore, the MOSFUN errors appear larger, relative to GENFIT, than they really are.

#### Potassium Ferrocyanide

The potassium ferrocyanide crystal is a 1 inch diameter by 1/8 inch thick sample (number E26532). A seven hour run was made to obtain data for analysis. The data was analyzed using Lorentzian lineshapes with the intensity model. Results show an isomer shift relative to natural iron of  $-0.0304$  mm/s compared to a reference value of  $-0.035 \pm 0.007$  mm/s (Ref 10: 154). The computed value is within the range of the standard

value. Chi-squared for this analysis was 2.045. The computed half-width was 0.127 mm/s. No other analysis was done with this sample.

#### Natural Iron

The natural iron absorber is a 25 micron thick foil provided with the cobalt-57 source. A four hour run was made to collect data and the data was analyzed using GENFIT. The CALFUN subroutine is called EAFECAL. The data was then evaluated using MOSFUN with the peak area as the intensity model (theory 2-Lorentzian and theory 4-Voigt). The areas of peaks 1, 2, and 3 were required to equal the areas of peaks 6, 5, and 4 respectively. Comparative results are shown in Table I. Figure 2 shows a representative plot of the fit obtained for peaks three and four (the inner two iron peaks) using the Lorentzian area as the intensity.

The results show a difference of less than 1% for all parameters except half-width. MOSFUN had a half-width that is 10% larger than GENFIT. The difference in the results obtained using the Voigt and Lorentzian lineshapes is negligible. However, the Voigt lineshape required 33 gradient iterations to converge compared to 7 for the Lorentzian.

All parameter values for either program are within 1% of the published standard values (Ref 10:154) for natural iron.

#### Sodium Nitroprusside

The sodium nitroprusside absorber is a National Bureau of Standards crystal measuring 1 cm by 1 cm by 0.0775 cm thick. A six hour spectrum was obtained and analyzed using

GENFIT. The CALFUN routine was NANPCAL. Analysis with MOSFUN was done using the peak intensity as the intensity parameter (theory 1 - Lorentzian and theory 3 - Voigt). Both peaks were required to have the same intensity and half-width. Results are presented in Table II, which compares the MOSFUN theories to GENFIT. Figure 3 shows a plot of the fit obtained using the Voigt lineshape.

The results show virtually no difference between GENFIT and either MOSFUN theory. Again the Voigt lineshape took longer than the Lorentzian to fit (41 gradient iterations to 22). Results for all methods are within 1% of the standard values (Ref 10:154).

#### Neodymium Cobalt

A neodymium cobalt-5 spectrum had been obtained in 1976 and was analyzed since it is a more complicated spectrum than any of the other samples used. Neodymium cobalt was the subject of Skluzacek's thesis. This particular spectrum was taken two years later after some equipment changes. The preparation of the sample is described by Skluzacek (Ref 8: 12-13) and will not be repeated.

Analysis with MOSFUN was done using theory model 1, Lorentzian intensities, and modeled three of the four sites modeled by Skluzacek (Ref 8:40-41). The fourth site has only a minor impact. In addition, the source used to obtain this spectrum was very weak and this reduced the relative intensities of the peaks by a factor of two, causing the fourth site to be buried. In the analysis, each site was assumed to

produce a six line iron spectrum with the intensities of peaks 1 and 3 held at a 3 to 1 ratio. The intensity of peak 2 was free to adjust.

The results were inconclusive due to the very weak intensity of the inner lines. The weak source presented poor statistics for the inner peaks of the stronger two sites and all of the third site. The outer peaks for the two stronger sites were converged upon and gave results that agreed well with Skluzacek's results. The stronger sites had magnetic fields of 239.9 kOe (kilo Oersted) and 204.7 kOe compared to 239.3 kOe and 205.5 kOe respectively for Skluzacek (Ref 8: 38-39). All other parameters are tied to the location of the inner peaks and no unique values could be found.

#### Program Considerations

MOSFUN is a user-oriented program and as such is very easy to use. However, parameter selection and control will make a great deal of difference in the results obtained. A TEST should be made to insure the initial selection of parameters and the correlation matrix are reasonable. The number of initial iterations should be kept small (approximately 5) and the parameters examined closely. The use of parameter limits will also keep parameters from deviating too far from reasonable values. To keep the number of free parameters low, use the same parameter for more than one peak whenever possible.

Table I  
Natural Iron Comparative Results : MOSFUN vs. GENFIT

Parameter	MOSFUN(Loren.)	MOSFUN(Voigt)	GENFIT
Chi-squared	2.15	2.35	2.33
Isomer Shift (mm/s)	-0.1054 $\pm$ .0024	-0.1055 $\pm$ .0024	-0.1057 $\pm$ .0005
Mag Field (kOe)	328.6 $\pm$ .4	328.5 $\pm$ .4	329.4 $\pm$ .1
Intensity Ratio			
peak 2 to 1	0.873 $\pm$ .014	0.875 $\pm$ .035	0.8708 $\pm$ .0098
peak 3 to 1	0.547 $\pm$ .086	0.549 $\pm$ .022	0.5451 $\pm$ .0088
Peak Separation(mm/s)			
peaks 1-6	10.571	10.571	10.578
peaks 2-5	6.121	6.123	6.125
peaks 3-4	1.671	1.669	1.672
Half-width (mm/s)			
peak 3	0.125 $\pm$ .006	0.124 $\pm$ .011	0.113 $\pm$ .005
peak 4	0.131 $\pm$ .007	0.131 $\pm$ .011	0.118 $\pm$ .011




Figure 2. Natural Iron Results for Peaks 3 and 4  
 MOSFET Using Lorentzian (Areas)

Table II  
Sodium Nitroprusside Comparative Results : MOSFUN vs. GENFIT

Parameter	MOSFUN(Loren.)	MOSFUN(Voigt)	GENFIT
Chi-squared	1.502	1.349	1.111
Isomer Shift (mm/s)	-0.367 $\pm$ 0.029	-0.367 $\pm$ 0.046	-0.3665 $\pm$ 0.0006
Quadrupole Split(mm/s)	1.692 $\pm$ 0.115	1.692 $\pm$ 0.181	1.693 $\pm$ 0.001
Peak Location (mm/s)			
peak 1	-1.213	-1.213	-1.213
peak 2	0.479	0.479	0.480
Half-width (mm/s)	0.139 $\pm$ 0.010	0.139 $\pm$ 0.018	0.137 $\pm$ 0.002




Figure 3. Sodium Nitroprusside Results  
 MOSFUN Using Voigt (Inten.)



## V. Conclusions and Recommendations

### Conclusions

The MOSFUN program greatly improves the analysis of Mossbauer spectra, even without the new hardware that is to be developed. MOSFUN is faster for the user, easier to control, more efficient due to this extra control, easier to expand, and gives results that agree within 1% with those obtained using GENFIT and a mainframe computer. MOSFUN will be able to analyze the data collected by the new hardware, when that hardware has been completed.

### Recommendations

The ease with which MOSFUN can be expanded and the parameters adjusted make this program the best to use for all future Mossbauer spectra evaluation. MOSFUN should be expanded to use the triangular wave mode and evaluated when triangular wave mode spectra can be generated. MOSFUN should be used to analyze spectra from the current MCA based system, even before the new hardware is developed. The data must be transferred to the disk manually, but this is still faster than the current evaluation procedure.

## Bibliography

1. Beller, D.E. "Crystallization Kinetics of Two Metallic Glasses by Mossbauer Spectroscopy." Unpublished MS thesis. Wright-Patterson AFB, Ohio: School of Engineering, Air Force Institute of Technology, March 1981.
2. Frauenfelder, Hans. The Mossbauer Effect. New York: W.A. Benjamin, Inc, 1962.
3. Muller, E.W. MOSFUN-A New and Versatile Mossbauer Fitting Program. Description and user's guide. Mainz, West Germany: Institut fur Anorganische Chemie und Analytische Chemie, Johannes Gutenberg-Universitat, 1980.  
[This document is available from the Mossbauer Effect Data Center]
4. Muller, E.W. "MOSFUN-A New and Versatile Mossbauer Fitting Program," Mossbauer Effect Reference and Data Journal, 4 (4): 89-90 (1 April 1981).
5. Roberts, L.D. "Mossbauer Studies of Metallic Glasses." Unpublished MS thesis. Wright-Patterson AFB, Ohio: School of Engineering, Air Force Institute of Technology, December 1978.
6. RT-11 System User's Guide. User's guide, Maynard, Massachusetts: Digital Equipment Corporation, March 1981. (AA-5279B-TC)
7. Schmidt, T.A. "An Analysis of Metallic Glasses by Mossbauer Spectroscopy." Unpublished MS thesis. Wright-Patterson AFB, Ohio: School of Engineering, Air Force Institute of Technology, March 1978.
8. Skluzacek, E.W. "Analysis of the Mossbauer Spectra of NdCo<sub>5</sub>." AFML-TR-75-162. Wright-Patterson AFB, Ohio: Air Force Materials Laboratory, 1976.
9. Skluzacek, E.W. Thesis Laboratory Notebook, Physics Department, Wright-Patterson AFB, Ohio, 1976.
10. Stevens, J.G. and Virginia E. Stevens. Mossbauer Effect Data Index: Covering the 1972 Literature. New York: IFI/PLENUM Data Corporation, 1973.
11. Towery, C.D. "Mossbauer Study of the Magnetocrystalline Anisotropy of NdCo<sub>5</sub>." Unpublished MS thesis. Wright-Patterson AFB, Ohio: School of Engineering, Air Force Institute of Technology, June 1971.

## Appendix A

### MOSFUN Installation Procedure

FORTRAN programs on the PDP 11 require three actions before they can be executed. Each FORTRAN routine is written and stored on the disk with a FOR file type. Each routine is compiled using the FORTRAN command (Ref 6:4-93) and an object file with an .OBJ file type created on the disk. Then the object files are linked together using the LINK command (Ref 6:4-119). LINK also adds any library routines that are needed. The final executable program is created with a .SAV file type. This is the only file needed for execution and is executed by a "RUN filename" command.

Three disk volumes are used to install and use MOSFUN. Volume 1 is the FORTRAN Working Volume and contains the FORTRAN compiler, operating system, and system library. It also contains the necessary utility programs to edit, compile, and link FORTRAN routines. Volume 2 is the MOSFUN Files Volume and contains all of the .FOR and .OBJ files to be linked to create MOSFUN.SAV. This volume is only used to edit and recompile MOSFUN if changes are made. Volume 3 is the MOSFUN Working Volume and contains MOSFUN.SAV and the data files necessary to analyze Mossbauer spectra using MOSFUN. This is the only volume needed to analyze spectra. Volumes 1 and 3 are bootable into the RT 11 single job operating system. Volume 2 is strictly a storage volume and must be used with one of the other volumes.

To edit a MOSFUN routine, put Volume 1 into drive DY0

and type "BOOT DY0:". This sets up the default values to edit and store routines on a second volume. Put Volume 2 into the other drive. Give the command "EDIT filename.FOR" and the appropriate FORTRAN routine will be put into the KED editor. After the changes have been made, compile the routine using the command "FORTRAN/NOLINENUMBER filename". The NOLINENUMBER option is required for MOSFUN to fit in the computer's memory. If errors exist and the line number is needed, the routine may be compiled with line numbers in order to correct errors. However, recompile without line numbers before linking. If a listing is desired, the option /LIST can be added to the FORTRAN command. If added before the routine name, the listing is sent to the lineprinter. If added after the routine name, the listing is stored on the disk with a .LST file type.

After a routine has been edited and successfully compiled, all the object files must be relinked. During the link procedure, the overlay regions are set up. To make the link procedure easier, an indirect command file "LNKMOS.COM" was set up. This file executes the commands needed to link all of the object files into the correct overlay regions and create MOSFUN.SAV. To link MOSFUN, give the command "LINK@LNKMOS". One major problem will frequently appear when this is done. Volume 2 is filled nearly to capacity. If there are not 300 contiguous blocks on the disk, a "SAV DEVICE FULL" error is generated due to the size of MOSFUN.SAV. If this occurs, delete the old .SAV file using the command "DELETE MOSFUN.SAV". Delete any backup text files (.BAK file types) and any other files that are not needed. Then squeeze the

remaining files so that all open blocks are at the end of the disk by giving the command "SQUEEZE DY1:". When this is completed try the link procedure again.

Once MOSFUN.SAV is ready, put Volume 3 into DY0 and boot it so that all default storage is to the MOSFUN Working Volume. Then copy MOSFUN.SAV using the command "COPY DY1: MOSFUN.SAV DY0:". MOSFUN is now ready to run again. Volumes 1 and 2 can be stored as they are no longer needed.

## Appendix B

### Subroutine DATI

Subroutine DATI controls the input of the experimental data from the disk. This subroutine was modified to identify and sort out counts from the internal time-keeping oscillator and the Moire interferometer for velocity measurement which are multiplexed with the experimental data by the spectrometer.

The identification text for the data, the format for the data, and several constant values are read from the file. The format for these inputs is described in Appendix E. If a format is not specified, a format of (8F10.0) is used for the experimental data. The first constant value was changed from the maximum velocity to the number of time-base overflows (TOVR). The next four constants still describe the number of data overflows (NOVF), the number of data points (N), the period or the number of channels a full period spectrum would contain (PERI), and the drive mode number (NDRV). A sixth and seventh constant were added to enter the number of channels between the channels containing Moire counts (MPLEX) and the first channel containing these counts (NST). All of these constants are transferred in common blocks, except NST which is transferred in the subroutine statement.

The time-base counts are read from the first 13 channels, overflows added, and then averaged so that the variable TIME contains the average number of time-base oscillator counts. These channels and channel 14 are zeroed

so they will not be considered in the rest of the program.  
TIME is transferred to subroutine DRIVE in common block FLD.

The velocity data (Moire counts) is read from channels NST to PERI/2 in jumps of MPLEX. Thus only the first half period is used for velocity calibration. If the spectrum is only a half period spectrum, then the entire spectrum is used since PERI/2 should equal N. The velocity points are checked for overflows, overflows are added, and the left half of the velocity curve made negative. These counts are then transferred to the array VCON. VCON(1) contains the number of velocity points (maximum of 64). VCON(2) to VCON(VCON(1)) contain the velocity data in counts which can be related to channel number in subroutine DRIVE.

A listing of the changed portion of subroutine DATI is attached. DATI after label 250 is unchanged from the original version.

# Listing of Subroutine DATI

```

SUBROUTINE DATI (N1,N2,N3)

C
COMMON /UNIT/ NRD,NWR,NCR,LPR,LFD,NPL,NVS,MAXN,MAXM,MAXF,NV(121)
COMMON /DAT/ VMAX,PERI,Y(1024),ID(40),NDRV,N,NF,VEL(10)
COMMON /THEO/ F,P(32),DF(32),NTEO,M,MF,KTEO(20),KP(32)
COMMON /FLD/ P3,B,V,VCON(65),V1,V2,GEO,B0,C0,NFLD,MSH,MPLX,TIN
COMMON /TRA/ IN(81),OUT(16)

C
DIMENSION KFRM(8),IFRM(8),INST(4)
DATA IFRM/2H(1,2H0F,2H0.,2H0),2H ,2H ,2H ,2H /
DATA INST(1)/1HL/,INST(2)/1HZ/,INST(3)/1HH/,INST(4)/1H*/
DATA NY/1HY/

C
C
C
C
N1 = INPUT UNIT
N2 = FACTOR FOR SEARCHING DROPPED CHANNELS
N3.NE.0 ERROR RETURN

N9=N2
IF(N1.EQ.0) GOTO 285
READ(N1,5200,END=202)(ID(I),I=1,40)
READ(N1,5200,END=202)(KFRM(I),I=1,8)
READ(N1,5210,END=202)(IN(I),I=1,80)
GOTO 203

C
C
C
C
EOF

202 WRITE(NWR,6205) N1
N=0
N3=1
RETURN

203 CALL TRANS(NN)
TOVR=OUT(1)
NOVF=OUT(2)
N=OUT(3)
PERI=OUT(4)
NDRV=OUT(5)
MPLEX=OUT(6)
NST=OUT(7)

205 IF (N.EQ.0) N=400

C
C
C
C
MORE THAN MAXN POINTS

IF(N.LE.MAX) : GOTO 220
N=0
N3=1
WRITE(NWR,6200) MAXN
RETURN

220 IF(NDRV.EQ.0) NDRV=3
IF(PERI.EQ.0.) PERI=N*2.

```



```

        IF(NST.EQ.0) NST=19
        P3=PERI * .5 + .5
230      IF(KFRM(1).NE.IFRM(8)) GOTO 240
        DO 231 I=1,8
231      KFRM(I)=IFRM(I)
240      WRITE(NWR,6230) (ID(I),I=1,39)
        WRITE(NWR,6240) TOVR,NOVF,N,PERI,NDRV,KFRM
        DO 241 I=1,MAXN
241      Y(I)=0.
        READ(N1,KFRM,END=250) (Y(I),I=1,N)
C
C      CALCULATE AVERAGE TIME OSC. COUNTS
C
99      TIME=0.
        DO 100 I=1,13
        TIME=TIME+Y(I)+10000000.*TOVR
100     Y(I)=0.
        TIME=TIME/13.
        Y(14)=0.
C
C      EXTRACT VELOCITY DATA, ACCOUNT FOR
C      OVERFLOWS AND ZERO AFTER PUTTING IN VCON.
C
        IIII=0
        III=0
        NEND=PERI/2
        DO 101 I=19,NEND,MPLEX
        VJMP=Y(I+MPLEX)-Y(I)
        IF(VJMP.GE.500000.) III=III+1
101     CONTINUE
        J=1
        DO 102 I=19,NEND,MPLEX
        J=J+1
        IF(III.EQ.0) GOTO 105
        VJMP=Y(I+MPLEX)-Y(I)
        IF(J.GT.65) GOTO 102
        VCON(J)=Y(I)+FLOAT(III)*1000000.
        IF(VJMP.GE.500000.) III=III-1
        GOTO 102
105     VJMP=Y(I)-Y(I+MPLEX)
        IF(J.GT.65) GOTO 102
        VCON(J)=Y(I)+FLOAT(IIII)*1000000.
        IF(VJMP.GE.500000.) IIII=IIII+1
102     Y(I)=0.
        Y(399)=0.
        III=-1
        DO 106 I=3,J-1
        VJMP=ABS(VCON(I-1))-VCON(I)
        IF(ABS(VCON(I)-VCON(I+1)).GT.(2.*VJMP)) III=1
        VCON(I)=VCON(I)*FLOAT(III)
106     CONTINUE
        VCON(2)=-VCON(2)
        VCON(1)=J

```

NN=N1  
NK=N2  
N1=0  
N2=NST  
CALL DRIVE(N1,N2)  
N1=NN  
N2=NK  
GOTO 260

C  
C  
C

LESS THAN N DATA IN PARTITION N1

250 WRITE(NWR,6210) I,N

## Appendix C

### Subroutine DRIVE

Subroutine DRIVE calculates the velocity, geometry factor, and derivatives of the first three parameters (baseline, geometry, and fold point) for any channel. The original subroutine DRIVE used the maximum velocity and the period to determine the velocity contribution to these values. Subroutine DRIVE was changed to allow a linear least squares fit to be performed on the velocity data extracted by subroutine DATI. This velocity fit is then used to find the above values.

Two arrays are built to use during the velocity fit. Array X(I) contains the channel number corresponding to the velocity in VCON(I). Array XY(I) contains VCON(I) which has been converted from counts to millimeters per second (mm/s) by the formula:

$$XY(I) = 317.5 \frac{VCON(I)}{TIME} \quad (C1)$$

The range of array X is scaled to a range of +1 to -1. The sums of powers of X (XDP(1-3)) and the sums of XY times the powers of X (XYDP(1-2)) are calculated. The values of XDP are placed in a symmetric array A and the values of XYDP placed in vector B. Subroutine is called to solve the system of linear equations  $AZ=B$ . The results (Z) and the range of X are scaled back to the original values. The velocity calibration is stored in the vector VEL. The velocity (V) in any channel (X) can then be calculated by the formula:

$$V(\text{mm/s}) = \text{VEL}(1) * X + \text{VEL}(2) \quad (\text{C2})$$

A listing of subroutine DRIVE is attached. Only the flyback mode calibration is shown, since it was the only drive mode used during this study.

# Listing of Subroutine DRIVE

```

SUBROUTINE DRIVE (N1,N2)
COMMON /UNIT/ NRD,NWR,NCR,LPR,LFD,NPL,NVS,MAXN,MAXM,MAXF,NV(12)
COMMON /DAT/ VMAX,PERI,Y(1024),ID(40),NDRV,N,NF,VEL(10)
COMMON /THEO/ F,P(32),DF(32),NTEO,M,MF,KTEO(20),KP(32)
COMMON /FLD/ P3,B,V,VCON(65),V1,V2,GEO,B0,C0,NFLD,MSH,MPLX,T
COMMON /TRA/ IN(81),OUT(16)
COMMON /LIN/ BB(32),A(32,32)

C
C   DIMENSION X(65),XY(65),XDP(10),XYDP(10)
C
C   N1 = 0  CALCULATE VELOCITY COEFF.
C   N1 = 1  VELOCITY V FROM CHANNEL NUMBER XX
C   N1 = 2  V,GEO
C   N1 = 3  DERIVATIVES P(1)-P(3)
C   N1 = 4  DERIVATIVES P(1)-P(2)
C
C   N2 = 1st CHANNEL CONTAINING VELOCITY (N1=0)
C   N2 = CHANNEL NUMBER   (FOR N1 = 1-4)
C   N2 = 0  ERROR RETURN
C
C   NDRV = 1  TRIANGULAR WAVE MODE
C   NDRV = 2  SINE WAVE MODE
C   NDRV = 3  FLYBACK WAVE MODE
C
C   IF(N1.NE.0) GOTO 1000
C   NST=N2
C   NFN=VCON(1)
C   GOTO (1,2,3), NDRV
1  WRITE(NWR,6000)
   RETURN
2  WRITE(NWR,6200)
   VCON(1)=0.
   RETURN
3  WRITE(NWR,6300)
   X(2)=NST
   XY(2)=VCON(2)*156.25/TIME
   DO 401 I=3,NFN
   X(I)=X(I-1)+MPLX
401 XY(I)=VCON(I)*156.25/TIME
   L=2
   K=1
C
C   FIND MIN AND MAX VALUES FOR X
C
C   XMIN=X(2)
C   XMAX=X(2)
C   DO 402 I=3,NFN
C   XMIN=AMIN1(XMIN,X(I))
402 XMAX=AMAX1(XMAX,X(I))

```

```

C
C      ZERO ARRAYS FOR SUMMING
C
      MM=2*K+1
      DO 403 I=1,MM
      XDP(I)=0.
      XYDP(I)=0.
403    CONTINUE
C
C      TRANSFORM RANGE OF X TO (-1,1) AND COMPUTE SUMS OF
C      POWERS OF X AND SUMS OF XY TIMES POWERS OF X.
C
      C1=2.0/(XMAX-XMIN)
      C2=(XMAX+XMIN)/(XMAX-XMIN)
      LU=2*K+1
      LL=K+2
      DO 404 II=2,NFN
      XP=1.
      XI=C1*X(II)-C2
      DO 405 I=1,L
      XDP(I)=XDP(I)+XP
      XYDP(I)=XYDP(I)+XP*XY(II)
405    XP=XP*XI
      DO 404 I=LL,LU
      XDP(I)=XDP(I)+XP
404    XP=XP*XI
      DO 406 I=1,L
      BB(I)=XYDP(I)
      DO 406 J=1,L
      A(I,J)=XDP(I+J-1)
406    CONTINUE
C
C      CALL LINEQ AND FIT VELOCITY DATA
C
      N1=0
      N2=2
      CALL LINEQ(N1,N2)
C
C      MOVE VELOCITY CALIB. TO VEL ARRAY
C
      DO 407 I=1,L
      JK=K-I+2
407    VEL(JK)=BB(I)
C
C      ADJUST COEFF. TO ORIGINAL RANGE OF X
C
      DO 408 I=1,K
      DO 408 J=1,I
408    VEL(J)=VEL(J)*C1
      C1=(XMAX+XMIN)/2.0
      DO 409 I=1,K
      MM=L-I+1
      DO 409 J=2,MM

```

```

409 VEL(J)=-C1*VEL(J-1)+VEL(J)
C
C COMPUTE MAXIMUM AND ROOT MEAN SQUARE
C ERRORS AND OUTPUT ERROR ANALYSIS
C
WRITE(NWR,6400)
EMAX=0.
SUM=0.
VMAX=0.
DO 410 I=2,NFN
YC=VEL(1)
DO 411 J=1,K
411 YC=YC*X(I)+VEL(J+1)
VMAX=AMAX1(VMAX,ABS(YC))
DIFF=YC-XY(I)
IF (I-1.GT.L) GOTO 413
WRITE(NWR,6410) I,X(I),XY(I),YC,DIFF,VEL(I-1)
GOTO 412
413 WRITE(NWR,6420) I,X(I),XY(I),YC,DIFF
412 EMAX=AMAX1(EMAX,ABS(DIFF))
410 SUM=SUM+DIFF**2
ERMS=SQRT(SUM/FLOAT(NFN-1))
WRITE(NWR,6430) EMAX,ERMS
RETURN
C
1000 GOTO (10,100,200),NDRV
C
C FLYBACK MODE
C
200 IF(N1.GE.3) GOTO 250
XX=N2
PER=PERI/4.
IF(XX.LT.PERI/2.) GOTO 232
XX=XX-PERI
C1=P3-PER-XX
S=1.
GOTO 240
232 C1=P3+PER-XX
S=-1.
240 C4=S*VMAX/PER
V=VEL(1)*XX+VEL(2)
IF(N1.EQ.1) RETURN
CX=V/VMAX
C2=S*(CX*CX-1.)
C5=1.+P(2)*C2
GEO=1./C5/C5
RETURN
C
C DERIVATIVES
C
250 MM=KP(1)
IF(MM.GT.0) DF(1)=GEO
MM=KP(2)

```

```

IF (MM.GT.0) DF (MM) = -2. *GEO*B0*C2/C5
IF (N1.EQ.4) RETURN
MM=KP(3)
IF (MM.GT.0) DF (MM) = GEO*((V/XX)*C0-B0*P(2)*V*4./C5/VMAX/PER)
RETURN

```

C

```

5000 FORMAT(80A1)
6300 FORMAT(1H ,4X,'DRIVE MODE FLYBACK WAVE')
6400 FORMAT(1H ,//,10X,36HPOLYNOMIAL LEAST SQUARE FIT
1ANALYSIS,/,4H I,6X,7HX-GIVEN,7X,7HY-GIVEN,6X,
28HY-FITTED,8X,5HERROR,10X,6HVEL(I),/)
6410 FORMAT(1H ,I3,4X,4(1PE10.3,4X),1PE13.6)
6420 FORMAT(1H ,I3,4X,4(1PE10.3,4X))
6430 FORMAT(1H ,9X,5HEMAX=,1PE15.6,9X,5HERMS=,1PE15.6)
END

```



## Appendix D

### Subroutine PLOT

Subroutine PLOT performs an alphanumeric plot on the terminal or the lineprinter. The plot will always contain the original data, but will also include the theoretical spectrum, if it has been calculated. The plot is performed over a range of channels specified by two parameters, V1 and V2. If these are zero, the default range will plot between channel 20 and the last channel (N).

The Y-axis represents intensity scaled relative to the baseline, so that the most intense dip within the channel range is full scale. The X-axis represents velocity in mm/s.

The plot is done as a one line text vector. Integer variables are used to calculate the relative intensity of the data point or the theory value, and scale them between 1 and 50. A symbol (+ for original data, \* for theory, and # for both) is placed in the text vector at the corresponding value in the vector, and all other points are blank. This text vector is then printed out. Every tenth channel also has the velocity value printed on the X-axis.

A listing of subroutine PLOT is attached.

# Listing of Subroutine PLOT

```

SUBROUTINE PLOT (N1,N2)
COMMON /UNIT/ NRD,NWR,NCR,LPR,LFD,NPL,NVS,MAXN,MAXM,MAXF,NV(12)
COMMON /DAT/ VMAX,PERI,Y(1024),ID(40),NDRV,N,NF,VEL(10)
COMMON /THEO/ F,P(32),DF(32),NTEO,M,MF,KTEO(20),KP(32)
COMMON /FDAT/ G,GP,GOPT,GOLD,GAM,YY(1024)
COMMON /FLD/ P3,B,V,VCON(65),V1,V2,GEO,B0,C0,NFLD,MSH,MPLX,T
COMMON /TRA/ IN(81),OUT(16)

*****
N1 = 6  LPR PLOT
N1 = 7  CONSOLE PLOT

N2 = 0  YY (THEORY and ORIG DATA)
N2 = 1  Y  (ORIG DATA)
*****

N3,N4  CONSOLE PLOT FROM CHANNEL N3 TO CHANNEL N4

DIMENSION NTXT(58)
DATA NBL/1H /,NSTAR/1H*/,NPLUS/1H+/,NPND/1H#/
K1=N1
K2=N2
N3=20
N4=N
IF(V1.NE.0.)N3=V1
IF(V2.NE.0.)N4=V2
IF(N3.LE.0) N3=20
IF(N4.GT.N) N4=N
WRITE(K1,6705)ID
IF(KTEO(1).NE.NBL) WRITE(K1,6706)KTEO
WRITE(K1,6707)N3,N4
WRITE(K1,6708) (ID(J),J=1,6)

SCALE Y AXIS

BASE=Y(N3)
IF(P(1).NE.0.)BASE=P(1)
YMAX=0.
YMIN=1.E10
DO 1200 I=N3,N4
G1=Y(I)
IF(N2.EQ.0) G1=YY(I)
IF(G1.EQ.0.) GOTO 1200
YMAX=AMAX1(YMAX,G1)
YMIN=AMIN1(YMIN,G1)
1200 CONTINUE
G2=YMAX-YMIN
IF(BASE.LE.G2) BASE=YMAX
WRITE(K1,6710) G2/BASE,.5*G2/BASE

```

```

C
C
C
      BEGIN PLOT

      N1=1
      N2=N3
      CALL DRIVE(N1,N2)
      WRITE(K1,6720) V
      K3=N3+1
      K4=N4-1
      DO 1250 I=K3,K4
      G1=Y(I)
      G6=YY(I)
      NQ=INT((G1-YMIN)/G2*50)
      MQ=INT((G6-YMIN)/G2*50)
      IF(G1.LE.0.) NQ=0
      IF(G6.LE.0.) MQ=0
      DO 1210 J=1,58
      NTXT(J)=NBL
      IF(J.EQ.NQ) NTXT(J)=NPLUS
      IF(J.EQ.MQ) NTXT(J)=NSTAR
      IF(J.EQ.MQ.AND.J.EQ.NQ) NTXT(J)=NPND
1210  CONTINUE
      IF(MOD(I,10).NE.0) GOTO 1220
      N1=1
      N2=I
      CALL DRIVE(N1,N2)
      WRITE(K1,6730) V,NTXT
      GOTO 1250
1220  WRITE(K1,6740) NTXT
1250  CONTINUE
      N1=1
      N2=N4
      CALL DRIVE(N1,N2)
      WRITE(K1,6720) V
      RETURN

C
6700  FORMAT(1H ,20X,'MOSFUN PLOT OF ',15X,'LEGEND FOR PLOT',/,22X,6I
      118X,'+ = ORIG DATA',/, ' VEL.',48X,'* = THEORY',/, '(MM/S)',48X,
      2' # = BOTH',/, ' V',20X,'INTENSITY RELATIVE TO BASELINE')
6705  FORMAT(1H ,//,5X,40A2)
6706  FORMAT(1H , 'THEORY= ',10A2,5X,10A2)
6707  FORMAT(1H ,/,15X,'CHANNEL ',I4,' TO CHANNEL ',I4)
6710  FORMAT(1H , ' V',4X,F7.3,19X,F7.3,22X,1HI)
6720  FORMAT(1H ,F7.3,'---I-----I
      1-----I-----I')
6730  FORMAT(1H ,F7.3,3X,58A1,1HI)
6740  FORMAT(1H ,5X,1HI,4X,58A1,1HI)
      END

```

## Appendix E

### MOSFUN User Supplement

This appendix contains supplemental information on how to enter the initial values required for MOSFUN to operate. This information adds to the MOSFUN documentation and is necessary due to the changes made to DRIVE. It also expands the information on the entry of the parameters and the correlation matrix.

MOSFUN requires three sets of input before a spectrum can be evaluated. These are the experimental data, the parameter set, and the correlation matrix. Examples of these inputs are attached. Example 1 is an experimental data disk file for sodium nitroprusside (SNP), example 2 is the interactive sequence to enter the parameter file for SNP, and example 3 is the interactive procedure to enter the correlation matrix for SNP. In examples 2 and 3, the computer printed entries are preceded by the identifier "CMPTR:". All other entries would be typed in by the user. These examples may be referred to during the procedure descriptions to follow.

The experimental data is entered from the disk along with several other values. These entries are an identification text for the data, the format by which the data is to be read, and constants for the number of time-base overflows (TOVR), the number of data overflows (NOVF), the total number of channels to be read (N), the number of channels in a full period spectrum (PERI), the number representing the drive

mode (NDRV), the multiplex interval for the velocity data (MPLEX), and the first channel containing velocity data.

The parameter set is entered from the disk or the console. Each parameter is composed of an identifying number, a name, a value, and upper and lower limits. Also entered with the parameters are the theory number and the total number of parameters.

The correlation matrix is entered from the disk or the console. The correlation matrix relates the parameters to the lines that make up the spectrum. Entries made with the correlation matrix are descriptive text to identify the correlation matrix, the total number of lines in the spectrum, the number of parameters that will be used to determine the intensity, position, half-width, and form of the lines, and how much of each parameter is to be applied to each line.

### Data

The data inputs are stored on the disk. They are read by MOSFUN when the command "DATIN" is given. MOSFUN asks for the location and file name for the data by the prompt "DYX:FILE.EXT". The file name (FILE) can have up to six characters but the file type (EXT) must have three characters. For example, to enter the sodium nitroprusside data from drive DY0, the response would be "DY0:SNP827.EXT".

A correlation factor will be requested by the prompt "CORR. FACTOR". The factor is entered from the console and is used to determine bad data points. A factor of 100 is sufficient for most sets of data.

The data is stored on the disk in any format, but the other information must precede it in three records. Each record can contain up to 80 characters and is ended by a carriage return. The first record contains text to identify the data. The second record contains the FORTRAN format that MOSFUN will use to read the data. The format must be inclosed in parentheses. The default format is (8F10.0) and is selected by a blank second record. The third record contains the constant values required to describe the data. The values are separated by commas and are entered in the following order:

1. TOVR - the number of time-base overflows
2. NOVF - the number of data overflows
3. N - the total number of channels to be read  
(default = 400)
4. PERI - the number of channels in a full period  
(default = N)
5. NDRV - the number representing the drive mode
6. MPLEX- the multiplex interval for the Moire counts
7. NST - the first channel containing Moire counts  
(default = 19)

Default values can be selected by entering a zero or leaving the field blank.

#### Parameter Set

The parameter set may be read from the disk or the console. Parameter entry is started by the command "PARIN" from the console. MOSFUN prompts "INPUT FROM FILE (Y/N)" to ask whether the input is coming from a disk file or the

terminal. If the parameters are to be read from a disk file, MOSFUN asks for the file location and name in the same way it did for the data. However, the parameter file name must be six characters long. If the parameters are to be entered manually or changes are to be made to the current parameter set, MOSFUN begins an interactive session (see example 2).

The example uses sodium nitroprusside and uses nine parameters. Besides the required three baseline parameters, the parameters are the intensity for each of the two sodium nitroprusside lines, the isomer shift, the quadrupole coupling, and the half-width of each line.

Parameters must be entered from the console initially. MOSFUN prompts with the current value, changes are made, and the changed values are restated. The process continues until stopped by the user. The first input is the theory number and the number of parameters. MOSFUN will prompt "THEORY 1 , 0 PARAMETERS". The theory number selected and the number of parameters are entered, separated by a comma. Each parameter's values are entered in a like manner. MOSFUN will prompt with the number and name, such as "NR 1, NAME BASE (\*=END)". A number and a name are required initially to define the parameter, but only the number is needed after it is defined. The program will then prompt with the current values for that parameter number in the form "VALUE 0.0000, LLIM 0, ULIM 0". The correct value and limits are entered, separated by commas. If no limits are entered, a value of zero is assigned and the parameter can vary without limits. If no values are entered, the previous values remain unchanged. All changed values are

restated by MOSFUN for review. The process continues until an astrisk is entered for a parameter number.

MOSFUN will then prompt "CHANGE SAVE FILE(Y/N)", to ask if the parameter set is to be saved to the disk. If not, PARIN stops. If yes, the program will ask which disk drive and what file name is to be used to store the file with the prompt "DYX:PAROUT.PAR". The response will state the drive, the file name, and the file type, and must contain 14 characters as in the prompt. The next prompt, "FILE NEW OR OLD", asks if the file is new or already on the disk. The file is then saved to the disk.

#### Correlation Matrix

The correlation matrix is read from the disk with the parameter set on the command "PARIN". To enter the correlation matrix from the console or to change any values, the command "CORR" is entered from the console, and an interactive procedure such as that in example 3 is started.

The first prompt asks where the final correlation matrix values are to be printed. This can be to the line-printer, the console, or not at all. This printout will be done after all values or changes have been made but before the changes are saved.

The first input to the correlation matrix file is a short (20 characters) text entry that describes the correlation matrix. MOSFUN prompts with "THEORY (20 CHAR)".

The next input to the file is five integers that tell the number of lines and the number of parameters that



describe the intensity, position, width, and form of the lines. MOSFUN prompts with the current values of these integers, and changes are made by entering the new values, separated by commas. All five must be entered, even if only one is changed. A blank input leaves all values unchanged.

The example used is sodium nitroprusside and has two of each type of parameter; intensity, position, and half-width. The intensity of each peak is allowed to vary independently as are the respective half-widths. The position of line 1 is the isomer shift less one-half of the quadrupole coupling. The position of line 2 is the isomer shift plus one-half of the quadrupole coupling. There are no form parameters since a Lorentzian lineshape is assumed.

The correlation matrix is entered next as four submatrices, one each for intensity, position, width, and form. Each has one row for each line in the spectrum and one column for each parameter that affects that submatrix. For example, if there are six lines and four intensity parameters, then the intensity submatrix will have six rows and four columns. MOSFUN tells which submatrix is to be entered by the prompt "PARAMETER TYPE I (I, P, G, X, \*=END)", where I is intensity, P is position, G is width, and X is form. Any submatrix can be selected by entering the appropriate letter. If the prompted entry is correct, a blank entry will begin the procedure to read the submatrix values. MOSFUN tells the column to be changed by a prompt that relates the submatrix column to a parameter. The prompt "1. COLUMN 4. PARAMETER NAME (\*=END)" means that column 1 relates parameter 4, called NAME, to the

lines in the spectrum. To change the column or parameter, enter the correct column number and the corresponding parameter number separated by a comma. A blank entry leaves the values unchanged. The values in that column are then listed, and changes can be made to all values in the column or to any row in the column. The values relate how many times that parameter will apply to each line. Once changes are made, an astrisk must be placed in the first column of an input to each of the above procedures.

Once all of the columns have been entered, the entire correlation matrix is printed as earlier directed. The prompt "CHANGE SAVE FILE (Y/N)" again appears and is exactly the same as for the parameter set, since both are stored in the same disk file.

# Example 1. SNP Data Input File

SNP(8-27-81) SODIUM NITROPRUSSIDE

5,0,0,400,000,3,10

557102	556975	557036	556968	557036	557058	556941	556996	557076	556989
556978	557065	556970	168526	113166	113203	113694	113592	316401	113550
113193	113216	113082	113254	113901	113079	113183	112830	298433	113523
113597	113546	113513	113255	113438	113667	113764	113178	281824	113116
113251	113060	113174	113549	112774	113385	113458	113172	265269	113320
112939	113135	113036	112953	113223	113861	113031	112738	247683	112383
113001	113101	112674	112217	113591	112437	113079	112768	230960	113213
112538	113052	112519	113105	113436	113215	113226	112919	213998	113104
113184	112865	113469	112378	113314	113364	112734	112787	196809	112938
112578	112983	112555	112611	112623	113663	112809	112845	179618	113079
112717	112506	112480	112809	113373	112785	112763	112544	162369	112951
112886	112229	112783	112609	112760	113093	112279	112923	144780	112760
112475	112940	112312	112953	112396	112598	112674	112594	127628	112437
112759	112203	112323	112920	113491	113003	113271	112707	110903	112397
112856	112814	112331	112261	112678	112820	112372	112809	93370	112778
112544	112458	112287	111393	112031	112281	111557	112316	76035	112388
112037	111446	111398	110862	111593	111450	111222	110199	59446	100923
108592	106309	105167	101910	97619	91400	84279	79871	42160	88437
95275	100570	103483	106098	107892	108989	110162	110054	25011	111205
110566	111043	110771	111137	111073	110620	111228	111021	8121	110724
110617	109949	110156	109316	109303	108100	106191	104924	10734	98517
91581	84673	80527	81037	88621	96133	100920	103413	25917	107300
109432	109410	110250	110901	110970	111012	110369	111398	45317	111436
111733	111776	111471	112079	111999	111912	111480	112200	61886	112133
112420	112310	112386	112906	112209	112012	112262	112687	78168	112764
112963	112509	112665	111829	112386	112674	112714	112786	96035	112148
112775	112381	112765	113000	112709	112654	112789	113040	112677	112510
112491	112748	112732	112827	112179	112861	112323	112546	130801	112958
112881	113030	112924	112508	112629	113006	112893	113054	147547	112328
113090	112301	112450	112840	112594	112189	112873	112564	164496	112588
113121	112711	113489	112620	113532	113002	113023	113366	182313	113070
113208	112889	112184	113348	112763	112688	112864	112882	199656	112855
112700	113055	112762	113085	112075	112548	113131	112866	216300	112997
112972	113144	113037	113043	112283	113204	113156	112907	233574	113024
113573	113488	113327	113108	113057	113189	113057	113613	250489	113625
113847	112638	113344	113281	113335	113273	113448	113351	267473	112535
113147	113153	113988	113522	112741	113354	113117	113624	284430	113328
112893	113944	112662	113041	113604	113197	113676	113718	300843	113165
113731	112910	113155	113096	113663	113717	112810	113106	318098	113895
113455	113724	113936	113627	113072	113853	113363	113132	334632	113930

# EXAMPLE 2 Entering Parameter File

```

CMPTR: *****COMMAND*****
      PA
CMPTR: *****PA*****
CMPTR: INPUT FROM FILE (Y/N)
      N

CMPTR: THEORY 1    0 PARAMETERS
      1,9

CMPTR: THEORY 1  LORENTZIAN (INTEN.)    9 PARAMETERS
CMPTR: NR 1    NAME BASE (*=END)
      hit carriage return (cr)

CMPTR: VALUE      0.00000,    LLIM    0.00000,    ULIM    0.00000
      123456

CMPTR:      1    BASE    123456.00000    0.00000    0.00000
CMPTR: NR 2    NAME GEO    (*=END)
      cr

CMPTR: VALUE      0.00000,    LLIM    0.00000,    ULIM    0.00000
      .001

CMPTR:      2    GEO      0.00100    0.00000    0.00000
CMPTR: NR 3    NAME FOLDP    (*=END)
      cr

CMPTR: VALUE      0.00000,    LLIM    0.00000,    ULIM    0.00000
      0

CMPTR:      3    FOLDP    0.00000    0.00000    0.00000
CMPTR: NR 4    NAME      (*=END)
      4 INTEN1

CMPTR: VALUE      0.00000,    LLIM    0.00000,    ULIM    0.00000
      32000

CMPTR:      4    INTEN1  32000.00000    0.00000    0.00000
CMPTR: NR 5    NAME      (*=END)
      5 INTEN2

CMPTR: VALUE      0.00000,    LLIM    0.00000,    ULIM    0.00000
      32000

CMPTR:      5    INTEN2  32000.00000    0.00000    0.00000
CMPTR: NR 6    NAME      (*=END)
      6 ISOSFT

CMPTR: VALUE      0.00000,    LLIM    0.00000,    ULIM    0.00000
      -.367

```

CMPTR:	6	ISOSFT	-0.36700	0.00000	0.00000
CMPTR:	NR 7	NAME	(*END)		
	7 QDSPL				
CMPTR:	VALUE	0.00000,	LLIM	0.00000,	ULIM 0.00000
	1.69				
CMPTR:	7	QDSPL	1.69000	0.00000	0.00000
CMPTR:	NR 8	NAME	(*END)		
	8 HFWD1				
CMPTR:	VALUE	0.00000,	LLIM	0.00000,	ULIM 0.00000
	.15, .1, .25				
CMPTR:	8	HFWD1	0.15000	0.10000	0.25000
CMPTR:	NR 9	NAME	(*END)		
	9 HFWD2				
CMPTR:	VALUE	0.00000,	LLIM	0.00000,	ULIM 0.00000
	.15, .1, .25				
CMPTR:	9	HFWD2	0.15000	0.10000	0.25000
CMPTR:	NR 9	NAME	HFWD2 (*END)		
	*				
CMPTR:	*****COMMAND*****				

### EXAMPLE 3 Entering Correlation Matrix

```

CMPTR: *****COMMAND*****
      CO
CMPTR: *****CO*****

CMPTR: CONTROL PRINTOUT 0=NO, 1=CONSOLE, 2=LINEPRINTER
      1

CMPTR: THEORY (20 CHAR.)
      INTEN,HW FREE

CMPTR: NO. OF LINES,INT.-,POS.-,WIDTH-,FURTHER PARAM.
CMPTR:      0      0      0      0      0
      2,2,2,2,0

CMPTR: PARAMETER TYPE I (I,P,G,X,*=END)
      hit carriage return (cr)

CMPTR:      1. COLUMN      4. PARAMETER = INTEN1 (*=END)
      cr

CMPTR:      1. COLUMN      4. PARAMETER = INTEN1 (*=END)
CMPTR:      1      0.00000
CMPTR:      2      0.00000
CMPTR: A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, *=NONE
      A

CMPTR: VALUES (MAX 16/LINE, *=END)
      1,0
      *

CMPTR:      1      1.00000
CMPTR:      2      0.00000
CMPTR: A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, *=END
      *

CMPTR:      2. COLUMN      5. PARAMETER = INTEN2 (*=END)
      cr

CMPTR:      2. COLUMN      5. PARAMETER = INTEN2 (*=END)
CMPTR:      1      0.00000
CMPTR:      2      0.00000
CMPTR: A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
      A

CMPTR: VALUES (MAX 16/LINE, *=END)
      0,1
      *

```

```

CMPTR:      1      0.00000
CMPTR:      2      1.00000
CMPTR:  A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
*

CMPTR:  PARAMETER TYPE P (I,P,G,X,*=END)
cr

CMPTR:      1. COLUMN      6. PARAMETER = ISOSFT (*=END)
cr

CMPTR:      1. COLUMN      6. PARAMETER = ISOSFT (*=END)
CMPTR:      1      0.00000
CMPTR:      2      0.00000
CMPTR:  A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
A

CMPTR:  VALUES (MAX 16/LINE, *=END)
1,1
*

CMPTR:      1      1.00000
CMPTR:      2      1.00000
CMPTR:  A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
*

CMPTR:      2. COLUMN      7. PARAMETER = QDSPL (*=END)
cr

CMPTR:      2. COLUMN      7. PARAMETER = QDSPL (*=END)
CMPTR:      1      0.00000
CMPTR:      2      0.00000
CMPTR:  A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
A

CMPTR:  VALUES (MAX 16/LINE, *=END)
-.5,.5
*

CMPTR:      1      -0.50000
CMPTR:      2      0.50000
CMPTR:  A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
*

CMPTR:  PARAMETER TYPE G (I,P,G,X,*=END)
cr

CMPTR:      1. COLUMN      8. PARAMETER = HFWD1 (*=END)
cr

CMPTR:      1. COLUMN      8. PARAMETER = HFWD1 (*=END)
CMPTR:      1      0.00000
CMPTR:      2      0.00000
CMPTR:  A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (*=END)
C

```

CMPTR: ROW NO., VALUE (\*=END)  
 1,1  
 2,0  
 \*

CMPTR: 1 1.000000  
 CMPTR: 2 0.000000  
 CMPTR: A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (\*=END)  
 \*

CMPTR: 2. COLUMN 9. PARAMETER = HFWD2 (\*=END)  
 cr

CMPTR: 2. COLUMN 9. PARAMETER = HFWD2 (\*=END)  
 CMPTR: 1 0.000000  
 CMPTR: 2 0.000000  
 CMPTR: A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (\*=END)  
 C

CMPTR: ROW NO., VALUE (\*=END)  
 2,1  
 \*

CMPTR: 1 0.000000  
 CMPTR: 2 1.000000  
 CMPTR: A=CHANGE ALL ELEMENTS, C=SINGLE ELEMENT, (\*=END)  
 \*

CMPTR: PARAMETER TYPE X (I,P,G,X,\*=END)  
 \*

CMPTR: THEORY INTEN,HW FREE  
 CMPTR:  
 CMPTR: 2 LINES  
 CMPTR:  
 CMPTR: 2 INTENSITY PARAMETERS  
 CMPTR: 4 5  
 CMPTR: 1 1.000000 0.000000  
 CMPTR: 2 0.000000 1.000000  
 CMPTR:  
 CMPTR: 2 POSITION PARAMETERS  
 CMPTR: 6 7  
 CMPTR: 1 1.000000 -0.500000  
 CMPTR: 2 1.000000 0.500000  
 CMPTR:  
 CMPTR: 2 HALF WIDTH PARAMETERS  
 CMPTR: 8 9  
 CMPTR: 1 1.000000 0.000000  
 CMPTR: 2 0.000000 1.000000  
 CMPTR:  
 CMPTR: 0 FURTHER PARAMETERS  
 CMPTR:  
 CMPTR: CHANGE SAVE FILE (Y/N)  
 Y



CMPTR: DYX:PAROUT.PAR (14 CHARACTERS)  
DYØ:SNP827.PAB

CMPTR: FILE NEW OR OLD  
NEW

CMPTR: \*\*\*\*\* COMMAND \*\*\*\*\*

### Vita

Bruce E. Pate was born on January 8, 1951 in Murfreesboro, Tennessee, and is the son of Gerald S. Pate and Donna J. Pate. He graduated from Quantico High School, Quantico, Virginia in June of 1969. In May of 1973, he graduated from Middle Tennessee State University, Murfreesboro, Tennessee with a Bachelor of Science Degree in Physics and Computer Science. He entered the Army through an ROTC commission in July of 1973. He has had duty at Fort Sheridan, Illinois and Aberdeen Proving Ground, Maryland as well as a tour in the Federal Republic of Germany. He and his wife, Kathy, have a daughter, Jenny, and a son, Mike.

Permanent Address: Route 7, Harrison Road  
Murfreesboro, Tennessee  
37130

Unclassified

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER AFIT/GNE/PH/82M-11	2. GOVT ACCESSION NO. AD-A124782	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Implementation of Software for a Microcomputer based Mossbauer Spectrometer		5. TYPE OF REPORT & PERIOD COVERED MS Thesis
7. AUTHOR(s) Bruce E. Pate		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS Air Force Institute of Technology (AFIT-EN) Wright-Patterson AFB, Ohio 45433		8. CONTRACT OR GRANT NUMBER(s)
11. CONTROLLING OFFICE NAME AND ADDRESS		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE March, 1982
		13. NUMBER OF PAGES 73
		15. SECURITY CLASS. (of this report) UNCLASSIFIED
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES Approved for public release; IAW AFR 190-17 LYNN E. WCLAVER Dean for Research and Professional Development Air Force Institute of Technology (AFIT) Wright-Patterson AFB, Ohio 45433 Director of Information 19 JAN 83		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Mossbauer Spectroscopy Computer Analysis Microcomputer programming		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) Computer analysis of Mossbauer spectra using a microcomputer requires tight controls on memory requirements and programming techniques. In this study, a FORTRAN IV program, MOSFUN, was programmed into an PDP-11 microcomputer. The program was modified to calculate the velocity for a data point using the least squares fit of counts from a Moire interferometer and to perform an alpha-numeric plot of the results and the original data. Spectra were		

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)

block 20 (cont)

obtained for absorbers of potassium ferrocyanide, natural iron, sodium nitroprusside, and neodymium cobalt. These spectra were then analyzed with the MOSFUN program and the results compared with analyses done with the former program, GENFIT, and with published values. Results showed MOSFUN to be more flexible, easier to adapt, and easier to control during the evaluation of spectra, while still providing results that agreed within 1% of those obtained with GENFIT and with standard values.

UNCLASSIFIED

END

FILMED

3-83

DTIC